



## REALQUALITY RQ-ADENO STANDARD

Ready-to-use quantified standards for quantification of *Adenovirus*,  
to be used on **GENEQUALITY** automatic platforms

### Code RQ-94-SA

Description	<p>The quantified REALQUALITY RQ-ADENO STANDARD is to be used exclusively in association with the REALQUALITY RQ-ADENO kit, code RQ-93, for quantitative testing on GENEQUALITY automatic platforms.</p> <p>The standard consists of DNA fragments containing the gene region of interest and therefore is not dangerous to the user.</p>
Reagent status	Ready-to-use
Packaging	<p>4 x 210 µL</p> <p>4 calibration points are provided, each sub-aliquoted in 3 x 70 µL, in sufficient quantity to do 9 standard curves.</p>
Stability	<p>12 months</p> <p>Each standard aliquote can undergo a maximum number of 3 freeze/thaw cycles: further cycles compromise the stability of the product and may alter its concentration.</p> <p>The standards cannot be used for more than 3 analytical runs.</p>
Storage	<p>-30°C/-20°C</p> <p>It is recommended to freeze the standard after each use.</p>
Concentration	<p>STANDARD 1: 20 viral genome copies /µL</p> <p>STANDARD 2: 200 viral genome copies /µL</p> <p>STANDARD 3: 2.000 viral genome copies /µL</p> <p>STANDARD 4: 20.000 viral genome copies /µL</p>
Instructions for use	<p>To obtain a standard curve in the range of 10<sup>2</sup>-10<sup>5</sup> copies of DNA, use 5 µL for reaction, strictly adhering to the indications given in the user manual of the REALQUALITY RQ-ADENO kit, code RQ-93.</p> <p>It is recommended to vortex and centrifuge each tube before use.</p>

### ORDER INFORMATION

Code	Product	PKG
RQ-94-SA	REALQUALITY RQ-ADENO STANDARD	4 x 210 µL
<b>In combination with the product:</b>		
RQ-93-4A	REALQUALITY RQ-ADENO	50 tests
RQ-93-6A		100 tests