



REALQUALITY STI-3 & STI-7

REALQUALITY STI-3 and STI-7 for reliable STI diagnostics across a broad range of clinically validated specimen types under IVDR.

Powered by GENEQUALITY® 2050, they provide a fully automated, standardized and traceable sample-to-result workflow.



**One Test. One Workflow.
One STI Solution.**

Real-Time
PCR



GOLD
STANDARD
DIAGNOSTICS



ANALITICA
A Gold Standard Diagnostics Company

STI-3
Single Multiplex PCR tube

STI-7
Two Multiplex PCR tubes



One workflow. One system.

One IVDR-compliant solution for STI diagnostics.

REALQUALITY
STI-7

STI-3

Chlamydia trachomatis (CT)

Neisseria gonorrhoeae (NG)

Mycoplasma genitalium (MG)

Trichomonas vaginalis (TV)

Mycoplasma hominis (MH)

Ureaplasma urealyticum (UU)

Ureaplasma parvum (UP)

**Dual Target
Region Amplification**

**Ready-to-use PCR mastermix
dUTP/UNG system**

Internal control (β-globin gene)

Positive controls included.

**One platform for comprehensive STI testing across
multiple specimen types**

Cervico-vaginal swab

Urethral swab

Urine




Seminal fluid

Rectal swabs


Oro-pharyngeal swabs

16/32

**Optimized throughput:
up to 32 STI-3 or 16
STI-7 samples per run**



**Vision system for
process and loading
monitoring**



**Magnetic bead
extraction with
integrated RT-PCR**

**from
150**

**Sample-to-result in
150 min for 16 samples,
165 min for 32 samples**

CODE	PRODUCT	FORMAT
RQ-196	REALQUALITY STI-7	108 tests on GENEQUALITY® platforms; 117 tests on validated thermal cyclers.
RQ-197	REALQUALITY STI-3	108 tests on GENEQUALITY® platforms; 117 tests on validated thermal cyclers.



AB ANALITICA srl: Via Svizzera, 16 - 35127 Padova ITALY | VAT No. 02375470289
Tel. +39 049 761689 | Fax. +39 049 8709510
www.abanalitica.com | customersupport@abanalitica.it

All promotional activities and marketing materials for AB ANALITICA devices are developed, reviewed, approved, and maintained in accordance with applicable regulatory requirements and supported by appropriate evidence. Marketing materials provide information that is consistent with the intended purpose, safety, and performance of the device and do not describe functions, features, or uses that are not supported by the applicable regulatory documentation. Complete and current information about the device is available in the latest version of the Instructions for Use (IFU). AB ANALITICA maintains documented procedures to ensure that promotional communications are accurate, comply with applicable regulatory requirements, and remain aligned with approved regulatory documentation throughout the device lifecycle. The information provided in this publication is subject to change without notice. For technical information and specific procedures, please contact your AB ANALITICA representative.